


Anti-GluR2/Gria2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P42262

Item Number: CM0006167



Introduction

Recombinant rabbit monoclonal antibody targeting human GluR2/GRIA2, validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience research on AMPA receptor-mediated synaptic transmission and plasticity.

[Learn More](#)

Parameter	Value
Product Name	GluR2/GRIA2 Rabbit mAb
Remarks/Alias	GLUR2; GLURB; GluA2; HBGR2; NEDLIB; glur-2; glur-B; GluR-K2; GluR2/GRIA2
Species	Human
GeneID	2891
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 330-470 of human GluR2/GRIA2 (NP_001077088.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42262
Protein Weight	99kDa
Shipping	Ice bag