

Anti-GluR4/GluA4/GRIA4 Rabbit Monoclonal Antibody For Wb And Elisa - P48058

Item Number: CM0012284

Introduction

Recombinant rabbit monoclonal antibody targeting human GluR4 (GRIA4, AMPA receptor subunit). Validated for Western blot and ELISA, cross-reacts with mouse and rat. Ideal for neuroscience research on excitatory synaptic transmission.

[Learn More](#)

Parameter	Value
Product Name	GluR4/GluA4/GRIA4 Rabbit mAb
Remarks/Alias	GLUR4; GLURD; GluA4; GLUR4C; NEDSGA; GluA4-ATD; GluR4/GluA4/GRIA4
Species	Human
Gene ID (Human)	2893
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P48058
Protein Weight	101kDa
Shipping	Ice bag