

Anti-GluR4/GluA4/GRIA4 Rabbit Polyclonal Antibody For Wb And Elisa - P48058

Item Number: CM0016010

Introduction

Rabbit polyclonal antibody targeting GluR4 (GRIA4), suitable for Western blotting and ELISA. Cross-reacts with human, mouse, and rat. Supports research on ionotropic glutamate receptors and fast excitatory synaptic transmission.

[Learn More](#)

Parameter	Value
Product Name	GluR4/GluA4/GRIA4 Rabbit pAb
Remarks/Alias	GLUR4; GLURD; GluA4; GLUR4C; NEDSGA; GluA4-ATD; GluR4/GluA4/GRIA4
Species	Human
GeneID (Human)	2893
GeneID	2893
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 254-433 of human GluR4/GluA4/GluR4/GluA4/GRIA4 (NP_001070712.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48058
Protein Weight	101kDa
Shipping	Ice bag