


Anti-Glut1/Slc2A1 Rabbit Monoclonal Antibody For If, Ihc-P, Elisa - P11166

Item Number: CM0003197



Introduction

Recombinant rabbit monoclonal antibody against human GLUT1/SLC2A1 (Glucose transporter type 1). Validated for IF, IHC-P, and ELISA with human, mouse, rat reactivity. Ideal for glucose metabolism, blood-brain barrier, and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	GLUT1/SLC2A1 Rabbit mAb
Remarks/Alias	CSE; PED; DYT9; GLUT; DYT17; DYT18; EIG12; GLUT1; HTLVR; GLUT-1; SDCHCN; GLUT1DS; GLUT1/SLC2A1
Species	Human
GeneID (Human)	6513
GeneID	6513
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11166
Protein Weight	54kDa
Shipping	Ice bag