


Anti-Glut1/Slc2A1 Rabbit Polyclonal Antibody For Wb And Elisa - P11166

Item Number: CM0014499



Introduction

Rabbit polyclonal antibody targeting human glucose transporter GLUT1/SLC2A1. Validated for western blot and ELISA, suitable for studying glucose transport, brain metabolism, and related diseases.

[Learn More](#)

Parameter	Value
Product Name	GLUT1/SLC2A1 Rabbit pAb
Remarks/Alias	CSE; PED; DYT9; GLUT; DYT17; DYT18; EIG12; GLUT1; HTLVR; GLUT-1; SDCHCN; GLUT1DS; GLUT1/SLC2A1
Species	Human
GeneID (human)	6513
GeneID	6513
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-400 of human GLUT1/SLC2A1 (NP_006507.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P11166
Protein Weight	54kDa
Shipping	Ice bag