


Anti-Glut1/Slc2A1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P11166

Item Number: CM0026171



Introduction

High-quality rabbit polyclonal antibody against GLUT1/SLC2A1, validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for glucose transport and metabolic disease research.

[Learn More](#)

Parameter	Details
Product Name	GLUT1/SLC2A1 Rabbit pAb
Remarks/Alias	CSE; PED; DYT9; GLUT; DYT17; DYT18; EIG12; GLUT1; HTLVR; GLUT-1; SDCHCN; GLUT1DS; GLUT1/SLC2A1
Species	Human
GeneID	6513
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 210-310 of human GLUT1/SLC2A1 (NP_006507.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11166
Protein Weight	54kDa
Shipping	Ice bag