

Anti-Glut2/Slc2A2 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P11168

Item Number: CM0031067

Introduction

Human GLUT2/SLC2A2 rabbit polyclonal antibody validated for WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for glucose transport and metabolic disease research in IVD settings.

[Learn More](#)

Parameter	Value
Product Name	GLUT2/SLC2A2 Rabbit pAb
Remarks/Alias	GLUT2; GLUT2/SLC2A2
Species	Human
Gene ID (Human)	6514
Gene ID	6514
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GLUT2/GLUT2/SLC2A2 (NP_000331.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11168
Protein Weight	57kDa
Shipping	Ice bag