


Anti-Glut3/Slc2A3 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P11169

Item Number: CM0008474



Introduction

High-quality anti-GLUT3/SLC2A3 rabbit monoclonal antibody validated for WB, IP, ELISA. Recognizes human, mouse, rat. Targets facilitated glucose transporter GLUT3 (P11169) with a molecular weight of 54 kDa.

[Learn More](#)

Parameter	Value
Product Name	GLUT3/SLC2A3 Rabbit mAb
Remarks/Alias	GLUT3; GLUT3/SLC2A3
Species	Human
Gene ID (Human)	6515
Gene ID	6515
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GLUT3/SLC2A3 (P11169).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11169
Protein Weight	54kDa
Shipping	Ice bag