


Anti-Glut4 Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P14672

Item Number: CM0008277



Introduction

GLUT4 Rabbit monoclonal antibody, validated for Western blot, immunofluorescence, IHC-P, and ELISA. Reacts with human, mouse, rat GLUT4. Ideal for diabetes and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	GLUT4 Rabbit mAb
Remarks/Alias	GLUT4
Species	Human
Gene ID (Human)	6517
Gene ID	6517
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 410-509 of human GLUT4 (NP_001033.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14672
Protein Weight	55 kDa/44 kDa
Shipping	Ice bag