


Anti-Glut4 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P14672

Item Number: CM0022135



Introduction

Rabbit polyclonal antibody targeting human GLUT4 (Glucose Transporter Type 4), suitable for WB, IF/ICC, IF-P, IHC-P, ELISA. Reacts with human and mouse. Ideal for glucose transport and insulin signaling studies.

[Learn More](#)

Parameter	Value
Product Name	GLUT4 Rabbit pAb
Remarks/Alias	GLUT4
Species	Human
GeneID	6517
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 210-310 of human GLUT4 (NP_001033.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P14672
Protein Weight	55kDa
Shipping	Ice bag