

Anti-Glutaminase (GLs) Rabbit Monoclonal Antibody For Wb - O94925

Item Number: CM0002283



Introduction

High-quality rabbit monoclonal antibody targeting human Glutaminase (GLS). Validated for WB, IF/ICC, and ELISA. Ideal for studying glutamine metabolism, kidney function, and neurotransmitter glutamate regulation. Immunogen: peptide aa 1-100 of human GLS.

[Learn More](#)

Parameter	Value
Product Name	Glutaminase (GLS) Rabbit mAb
Remarks/Alias	GAC; GAM; KGA; GLS1; AAD20; DEE71; GDPAG; CASGID; EIEE71; Glutaminase (GLS)
Species	Human
GeneID (human)	2744
GeneID	2744
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Glutaminase (GLS) (O94925)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O94925
Protein Weight	73kDa
Shipping	Ice bag