

Anti-Glutaminase (GLs) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O94925

Item Number: CM0028994



Introduction

Rabbit polyclonal antibody against human Glutaminase (GLS) validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for investigating glutamine metabolism and neurotransmitter regulation.

[Learn More](#)

Parameter	Value
Product Name	Glutaminase (GLS) Rabbit pAb
Remarks/Alias	GAC; GAM; KGA; GLS1; AAD20; DEE71; GDPAG; CASGID; EIEE71; Glutaminase (GLS)
Species	Human
GeneID (Human)	2744
GeneID	2744
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94925
Protein Weight	73kDa
Shipping	Ice bag