


Anti-Glutamine Synthetase (GluI) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P15104

Item Number: CM0006680



Introduction

High-quality rabbit monoclonal antibody against Glutamine Synthetase (GLUL), validated for WB, IHC-P, ELISA in human, mouse, rat. Ideal for glutamine metabolism and vascular development research.

[Learn More](#)

Parameter	Specification
Product Name	Glutamine Synthetase (GLUL) Rabbit mAb
Remarks/Alias	GS; GLNS; PIG43; PIG59; Glutamine Synthetase (GLUL)
Species	Human
GeneID (Human)	2752
GeneID	2752
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Glutamine Synthetase (GLUL) (P15104).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15104
Protein Weight	42kDa
Shipping	Ice bag