


Anti-Glutamine Synthetase (GluI) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P15104

Item Number: CM0022117



Introduction

High-quality rabbit polyclonal antibody against human glutamine synthetase (GLUL). Validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat GS. Ideal for nitrogen metabolism and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Glutamine Synthetase (GLUL) Rabbit pAb
Remarks/Alias	GS; GLNS; PIG43; PIG59; Glutamine Synthetase (GLUL)
Species	Human
GeneID (Human)	2752
GeneID	2752
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-373 of human Glutamine Synthetase (GLUL) (NP_002056.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15104
Protein Weight	42kDa
Shipping	Ice bag