

Anti-Glutathione Synthetase (Gss) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P48637

Item Number: CM0003537



Introduction

High-specificity rabbit monoclonal antibody against human Glutathione Synthetase (GSS), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for glutathione metabolism and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	Glutathione Synthetase (GSS) Rabbit mAb
Remarks/Alias	GSHS; HEL-S-64p; HEL-S-88n; Glutathione Synthetase (GSS)
Species	Human
Gene ID (Human)	2937
Gene ID	2937
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 350-450 of human Glutathione Synthetase (GSS) (P48637).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48637
Protein Weight	52kDa
Shipping	Ice bag