

## Anti-Glyat Polyclonal Antibody For Wb, Elisa - Q6Ib77

Item Number: CM0014456

### Introduction

GLYAT Rabbit pAb is a polyclonal antibody validated for WB and ELISA. It recognizes human GLYAT and cross-reacts with mouse and rat, making it suitable for detoxification, metabolism, and liver cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | GLYAT Rabbit pAb                                                                                                  |
| Remarks/Alias    | CAT; GAT; ACGNAT; GLYAT                                                                                           |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 10249                                                                                                             |
| Gene ID          | 10249                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-163 of human GLYAT (NP_005829.3). |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Mouse, Rat                                                                                                        |
| SWISS            | Q6IB77                                                                                                            |
| Protein Weight   | 34kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |