

Anti-Glyat/Gat Monoclonal Antibody For Wb, Ihc-P, Elisa - Q6Ib77

Item Number: CM0006889

Introduction

Rabbit monoclonal antibody targeting human GLYAT/GAT (Gene ID: 10249), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for metabolic detoxification and liver cancer research.

[Learn More](#)

Parameter	Value
Product Name	GLYAT/GAT Rabbit mAb
Remarks/Alias	CAT; GAT; ACGNAT; GLYAT/GAT
Species	Human
GeneID (human)	10249
GeneID	10249
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human GLYAT/GAT (NP_964011.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6IB77
Protein Weight	34kDa
Shipping	Ice bag