

Anti-Glycerol Kinase (Gk) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P32189

Item Number: CM0002676

Introduction

Glycerol kinase (GK) rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Targets human, mouse, rat GK (P32189), a key enzyme in glycerol metabolism and glycerolipid synthesis. Suitable for diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	Glycerol kinase (GK) Rabbit mAb
Remarks/Alias	GK1; GKD; Glycerol kinase (GK)
Species	Human
GeneID (Human)	2710
GeneID	2710
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32189
Protein Weight	61kDa
Shipping	Ice bag