

Anti-Glycerol Kinase (Gk) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P32189

Item Number: CM0027675

Introduction

Rabbit polyclonal antibody against human glycerol kinase (GK). Recognizes human, mouse, rat. Validated for WB, IHC-P, ELISA. Ideal for studies of glycerol metabolism, glycerolipid synthesis, and related metabolic pathways.

[Learn More](#)

Parameter	Value
Product Name	Glycerol kinase (GK) Rabbit pAb
Remarks/Alias	GK1; GKD; Glycerol kinase (GK)
Species	Human
GenelD (human)	2710
GenelD	2710
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 422-559 of Glycerol kinase (GK) (NP_001191948.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32189
Protein Weight	61kDa
Shipping	Ice bag