


Anti-Glycogen Synthase 1 (Gys1) Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P13807

Item Number: CM0021776



Introduction

Rabbit polyclonal antibody against human GYS1, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Useful for glycogen biosynthesis and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	Glycogen synthase 1 (GYS1) Rabbit pAb
Remarks/Alias	GSY; GYS; Glycogen synthase 1 (GYS1)
Species	Human
GeneID (human)	2997
GeneID	2997
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 488-737 of human Glycogen synthase 1 (NP_002094.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13807
Protein Weight	84kDa
Shipping	Ice bag