


Anti-Glycophorin A (Gypa) Rabbit Polyclonal Antibody For Wb And Elisa - P02724

Item Number: CM0021151



Introduction

Rabbit polyclonal antibody targeting human GYPa (Glycophorin A, CD235a, MN sialoglycoprotein). Validated for WB and ELISA. Major erythrocyte membrane protein, carries MN blood group antigens, and serves as receptor for HAV and P. falciparum. Molecular weight: ~16 kDa.

[Learn More](#)

Parameter	Details
Product Name	GYPa Rabbit pAb
Remarks/Alias	MN; GPA; MNS; GPSAT; PAS-2; CD235a; GP ^E rik; HGpMiv; HGpMIXI; HGpSta(C); GYPa
Species	Human
GeneID (Human)	2993
GeneID	2993
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-91 of human GYPa (P02724).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P02724
Protein Weight	16kDa
Shipping	Ice bag