


Anti-Glycophorin C (Gypc) Polyclonal Antibody For Wb, Elisa - P04921

Item Number: CM0017725



Introduction

Rabbit polyclonal antibody targeting Glycophorin C (GYPC), validated for Western blot and ELISA. Detects human and mouse GYPC (14 kDa), a sialoglycoprotein involved in erythrocyte stability and malaria receptor function.

[Learn More](#)

Parameter	Value
Product Name	Glycophorin C (GYPC) Rabbit pAb
Remarks/Alias	GE; GPC; GPD; GYPD; CD236; PAS-2; CD236R; PAS-2'; Glycophorin C (GYPC)
Species	Human
Gene ID (Human)	2995
Gene ID	2995
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-128 of human Glycophorin C (Glycophorin C (GYPC)) (NP_002092.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P04921
Protein Weight	14kDa
Shipping	Ice bag