

Anti-Glyctk Polyclonal Antibody For Wb, Elisa - Q8lvs8

Item Number: CM0024896

Introduction

High-quality rabbit polyclonal antibody against human GLYCTK (glycerate kinase), suitable for WB and ELISA. Reacts with human, mouse, and rat. Recombinant immunogen targeting C-terminal region. Ideal for studying mitochondrial metabolism and related research.

[Learn More](#)

Parameter	Value
Product Name	GLYCTK Rabbit pAb
Remarks/Alias	HBEBP2; HBEBP4; HBeAgBP4A; GLYCTK
Species	Human
Gene ID (Human)	132158
Gene ID	132158
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 284-523 of human GLYCTK (NP_660305.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IVS8
Protein Weight	55kDa
Shipping	Ice bag