

Anti-Glypican 3 (Gpc3) Monoclonal Antibody For Flow Cytometry - P51654

Item Number: CM0012580



Introduction

-conjugated rabbit monoclonal antibody against human Glypican 3 (GPC3), validated for flow cytometry. Specific to human GPC3, a cell surface proteoglycan involved in Wnt and Hedgehog signaling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Glypican 3 (GPC3) mAb
Remarks/Alias	SGB; DGSX; MXR7; SDYS; SGBS; OCI-5; SGBS1; GTR2-2
Species	Human
Gene ID (Human)	2719
Gene ID	2719
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 481-580 of human Glypican 3 (GPC3) (NP_004475.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P51654
Protein Weight	66kDa
Shipping	Ice bag