


Anti-Glypican 3 (Gpc3) Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - P51654

Item Number: CM0006802



Introduction

Rabbit monoclonal antibody against human Glypican 3 (GPC3), validated for WB, IF/ICC, FC, ELISA. Suitable for cancer, developmental biology, and Wnt/Hedgehog signaling research.

[Learn More](#)

Parameter	Value
Product Name	Glypican 3 (GPC3) Rabbit mAb
Remarks/Alias	SGB; DGSX; MXR7; SDYS; SGBS; OCI-5; SGBS1; GTR2-2; Glypican 3 (GPC3)
Species	Human
GeneID (human)	2719
GeneID	2719
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 481-580 of human Glypican 3 (GPC3) (NP_004475.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P51654
Protein Weight	66kDa
Shipping	Ice bag