


Anti-Glypican 3 (Gpc3) Polyclonal Antibody For Wb, Fc, Elisa - P51654

Item Number: CM0019017



Introduction

High-quality rabbit polyclonal antibody against Glypican 3 (GPC3) for Western blot, flow cytometry, and ELISA. Cross-reacts with human, mouse, and rat samples. Ideal for cancer and developmental biology research.

[Learn More](#)

Parameter	Specification
Product Name	Glypican 3 (GPC3) Rabbit pAb
Remarks/Alias	SGB; DGSX; MXR7; SDYS; SGBS; OCI-5; SGBS1; GTR2-2; Glypican 3 (GPC3)
Species	Human
GeneID (Human)	2719
GeneID	2719
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51654
Protein Weight	66kDa
Shipping	Ice bag