


Anti-Glyr1 Polyclonal Antibody For Wb, Elisa - Q49A26

Item Number: CM0025623



Introduction

Rabbit polyclonal antibody against GLYR1, validated for Western blot (WB) and ELISA in human, mouse, and rat. Ideal for chromatin biology, nucleosome remodeling, and gene expression regulation studies.

[Learn More](#)

Parameter	Value
Product Name	GLYR1 Rabbit pAb
Remarks/Alias	NP60; NPAC; hNDF; BM045; HIBDL; N-PAC; GLYR1
Species	Human
Gene ID	84656 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-270 of human GLYR1 (NP_115958.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q49A26
Protein Weight	61kDa
Shipping	Ice bag