

Anti-Glyt2 Monoclonal Antibody For Wb, Elisa - Q9Y345

Item Number: CM0002626

Introduction

Anti-Glyt2 Rabbit monoclonal antibody for WB, ELISA applications. Validated for Human, Mouse, Rat samples. Targets sodium- and chloride-dependent glycine transporter 2, critical for glycinergic neurotransmission. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Glyt2 Rabbit mAb
Remarks/Alias	NET1; GLYT2; HKPX3; GLYT-2; Glyt2
Species	Human
GeneID (Human)	9152
GeneID	9152
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 250-350 of human Glyt2 (Q9Y345)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y345
Protein Weight	87kDa
Shipping	Ice bag