

Anti-Gm-Csf Monoclonal Antibody For Wb, If/Icc, Elisa - P01587

Item Number: CM0011275

Introduction

Rabbit monoclonal antibody against mouse GM-CSF, validated for Western blot, immunofluorescence, and ELISA. Ideal for cytokine studies and hematopoietic differentiation research. SWISS: P01587.

[Learn More](#)

Parameter	Value
Product Name	GM-CSF Rabbit mAb
Remarks/Alias	CSF; Csfgm; GMCSF; Gm-CSf; MGI-IGM
Species	Mouse
Gene ID (Human)	1437
Gene ID (Mouse)	12981
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-141 of mouse GM-CSF (NP_034099.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	P01587
Protein Weight	16kDa
Shipping	Ice bag