

Anti-Gm-Csf Rabbit Monoclonal Antibody For Wb, Elisa - P48750

Item Number: CM0008688

Introduction

Recombinant rabbit monoclonal antibody targeting rat Granulocyte-macrophage colony-stimulating factor (GM-CSF, CSF2). Validated for WB and ELISA applications. 17 kDa secreted cytokine. For hematopoiesis research.

[Learn More](#)

Parameter	Value
Product Name	GM-CSF Rabbit mAb
Remarks/Alias	Gmcsf; Gm-csf
Species	Rat
GeneID (Human)	1437
GeneID	116630
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-144 of rat GM-CSF (NP_446304.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P48750
Protein Weight	17kDa
Shipping	Ice bag