

Anti-Gm130 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q08379

Item Number: CM0024600



Introduction

High-quality Anti-GM130 Rabbit Polyclonal Antibody (Q08379) validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat GM130 (Golgin A2). Ideal for studying Golgi structure, mitosis, and centrosome organization.

[Learn More](#)

Parameter	Value
Product Name	GM130 Rabbit pAb
Remarks/Alias	GM130; DEDHMB
Species	Human
GeneID (human)	2801
GeneID	2801
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human GM130 (NP_004477.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08379
Protein Weight	113kDa
Shipping	Ice bag