

Anti-Gm13125/Pramef20 Rabbit Polyclonal Antibody For Wb, IHC-P, Elisa - Q5Vt98

Item Number: CM0026842

Introduction

High-quality rabbit polyclonal antibody targeting mouse Gm13125 (PRAMEF20), suitable for WB, IHC-P, and ELISA applications. Exhibits cross-reactivity with human and rat, detecting a ~55 kDa band. Ideal for studying PRAME family proteins with leucine-rich repeats.

[Learn More](#)

Parameter	Value
Product Name	Gm13125 Rabbit pAb
Remarks/Alias	Gm13125; EG627009; Pramef20
Species	Mouse
GeneID (human)	645425
GeneID	627009
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 67-274 of mouse Gm13125 (NP_001108549.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5VT98
Protein Weight	55kDa
Shipping	Ice bag