

Anti-Gmfb Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P60983

Item Number: CM0028090

Introduction

Rabbit polyclonal antibody targeting human GMFB, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat, ideal for neural differentiation and cancer research.

[Learn More](#)

Parameter	Value
Product Name	GMFB Rabbit pAb
Remarks/Alias	GMF; GMFB
Species	Human
Gene ID (Human)	2764
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human GMFB (NP_004115.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60983
Protein Weight	17kDa
Shipping	Ice bag