

Anti-Gmip Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9P107

Item Number: CM0021256



Introduction

High-quality rabbit polyclonal antibody against human GMIP (GEM-interacting protein, ARHGAP46). Validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying RhoA GTPase regulation and cellular signaling.

[Learn More](#)

Parameter	Value
Product Name	GMIP Rabbit pAb
Remarks/Alias	ARHGAP46; GMIP
Species	Human
GeneID (human)	51291
GeneID	51291
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-91 of human GMIP (NP_057657.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P107
Protein Weight	107kDa
Shipping	Ice bag