


Anti-Gmn Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O75496

Item Number: CM0006596



Introduction

High-quality rabbit monoclonal antibody targeting human GMNN (Geminin). Validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human and mouse GMNN, a key inhibitor of DNA replication and cell cycle regulator.

[Learn More](#)

Parameter	Value
Product Name	GMNN Rabbit mAb
Remarks/Alias	Gem; MGORS6; Geminin
Species	Human
GeneID (human)	51053
GeneID	51053
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-209 of human GMNN (NP_056979.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75496
Protein Weight	24kDa
Shipping	Ice bag