

Anti-Gmps Polyclonal Antibody For Wb, Ip, Elisa - P49915

Item Number: CM0030795



Introduction

Rabbit polyclonal antibody targeting human GMPS (GMP synthase). Validated for WB, IP, ELISA; cross-reacts with mouse and rat. Immunogen: recombinant fusion protein corresponding to amino acids 394-693 of human GMPS.

[Learn More](#)

Parameter	Value
Product Name	GMPS Rabbit pAb
Remarks/Alias	GATD7; GMPS
Species	Human
GeneID (Human)	8833
GeneID	8833
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 394-693 of human GMPS (NP_003866.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49915
Protein Weight	77kDa
Shipping	Ice bag