

Anti-Gn1A Polyclonal Antibody For Wb, Elisa - Q4Adv8

Item Number: CM0022489

Introduction

Anti-Gn1a (OsCKX2) polyclonal antibody validated for WB and ELISA. Targets cytokinin dehydrogenase 2 in rice, a key regulator of grain yield. Sourced from rabbit; cross-reacts with *Oryza sativa*.

[Learn More](#)

Parameter	Value
Product Name	Gn1a Rabbit pAb
Remarks/Alias	CKX2; Gn1a; OsCKX2; Osj_00744; B1046G12.8
Species	Oryza sativa
GeneID	4327333
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-500 of Oryza sativa Gn1a (Q4ADV8)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q4ADV8
Protein Weight	60kDa
Shipping	Ice bag