

Anti-Gna11 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P29992

Item Number: CM0026792



Introduction

High-quality rabbit polyclonal antibody targeting GNA11 (G protein subunit alpha 11, SWISS P29992). Validated for WB, IHC-P, ELISA in human, mouse, rat. Ideal for GPCR signaling studies.

[Learn More](#)

Parameter	Value
Product Name	GNA11 Rabbit pAb
Remarks/Alias	FBH, FBH2, FHH2, HG1K, HHC2, GNA-11, HYPOC2
Species	Human
GeneID (Human)	2767
GeneID	2767
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-350 of human GNA11 (NP_002058.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29992
Protein Weight	42kDa
Shipping	Ice bag