


Anti-Gna13 Polyclonal Antibody For Wb, If/Icc, Elisa - Q14344

Item Number: CM0026448



Introduction

Rabbit polyclonal antibody targeting GNA13 (G protein subunit alpha-13). Validated for WB, IF/ICC, ELISA with human, mouse, rat cross-reactivity. Suitable for studying Rho signaling and cancer invasion mechanisms.

[Learn More](#)

Parameter	Description
Product Name	GNA13 Rabbit pAb
Remarks/Alias	G13; HG1N; GNA13
Species	Human
GeneID (Human)	10672
GeneID	10672
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human GNA13 (NP_006563.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14344
Protein Weight	44kDa
Shipping	Ice bag