


Anti-Gna13 Rabbit Polyclonal Antibody For Ihc-P - Q14344

Item Number: CM0019253



Introduction

Rabbit polyclonal antibody targeting GNA13 (G protein subunit alpha 13) validated for IHC-P in human, mouse, and rat tissues. Recombinant full-length immunogen; suitable for investigations of GPCR-RhoA signaling, cancer metastasis, and germinal center biology.

[Learn More](#)

Parameter	Value
Product Name	GNA13 Rabbit pAb
Remarks/Alias	GNA13; G13; G protein subunit alpha 13
Species	Human
Gene ID (Human)	10672
Gene ID	10672
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-377 of human GNA13 (NP_006563.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14344
Protein Weight	33kDa/44kDa
Shipping	Ice bag