

Anti-Gna13 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14344

Item Number: CM0021917



Introduction

Rabbit polyclonal antibody against human GNA13, validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, and rat GNA13. Ideal for studying G protein signaling, RhoA activation, and cancer metastasis.

[Learn More](#)

Parameter	Value
Product Name	GNA13 Rabbit pAb
Remarks/Alias	G13; HG1N; GNA13
Species	Human
GeneID (human)	10672
GeneID	10672
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human GNA13 (NP_006563.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14344
Protein Weight	44kDa
Shipping	Ice bag