

Anti-Gna15 Rabbit Polyclonal Antibody For Wb, Elisa - P30679

Item Number: CM0024358

Introduction

Rabbit polyclonal antibody targeting human GNA15 (G protein subunit alpha-15), suitable for WB and ELISA, with cross-reactivity to mouse and rat. Useful for studying GPCR signaling in hematopoietic cells.

[Learn More](#)

Parameter	Value
Product Name	GNA15 Rabbit pAb
Remarks/Alias	HG1L; GNA16; GNA15
Species	Human
Gene ID (Human)	2769
Gene ID	2769
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-374 of human GNA15 (NP_002059.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P30679
Protein Weight	44kDa
Shipping	Ice bag