

Anti-Gna15 Rabbit Polyclonal Antibody For Wb, Elisa - P30679

Item Number: CM0021404

Introduction

High-quality GNA15 rabbit polyclonal antibody validated for western blot and ELISA, targeting human and mouse G alpha-16 protein. Ideal for GPCR signaling, hematopoiesis, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	GNA15 Rabbit pAb
Remarks/Alias	HG1L; GNA16; GNA15
Species	Human
GeneID (human)	2769
GeneID	2769
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 112-374 of human GNA15 (NP_002059.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P30679
Protein Weight	44kDa
Shipping	Ice bag