


Anti-Gnai1 Rabbit Polyclonal Antibody For Wb, Elisa - P63096

Item Number: CM0027553



Introduction

Rabbit polyclonal antibody targeting GNAI1 (Gi, HG1B, NEDHISB). Detects human, mouse, rat. Suitable for Western blot and ELISA. Involved in GPCR signaling, inhibiting adenylate cyclase. Protein weight 40 kDa.

[Learn More](#)

Parameter	Value
Product Name	GNAI1 Rabbit pAb
Remarks/Alias	Gi; HG1B; NEDHISB; GNAI1
Species	Human
Gene ID (Human)	2770
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human GNAI1 (NP_002060.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63096
Protein Weight	40kDa
Shipping	Ice bag