


Anti-Gnai2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P04899

Item Number: CM0015191



Introduction

Rabbit polyclonal antibody against GNAI2 (G protein subunit alpha i2) validated for WB, IF/ICC, ELISA in human, mouse, rat. Suitable for studying GPCR signaling, angiogenesis, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	GNAI2 Rabbit pAb
Remarks/Alias	GIP; HG1C; GNAI2B; H_LUCA15.1; H_LUCA16.1; GNAI2
Species	Human
Gene ID (Human)	2771
Gene ID	2771
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-355 of human GNAI2 (NP_002061.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04899
Protein Weight	40kDa
Shipping	Ice bag