

Anti-Gnai3 Rabbit Polyclonal Antibody For Wb And Elisa - P08754

Item Number: CM0014935



Introduction

Rabbit polyclonal antibody against human GNAI3, validated for WB and ELISA. Reacts with human, mouse, and rat. Targets N-terminal epitope. For GPCR signaling and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	GNAI3 Rabbit pAb
Remarks/Alias	87U6, HG1A, ARCND1, ARCODS, GNAI3
Species	Human
Gene ID	2773
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human GNAI3 (NP_006487.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08754
Protein Weight	41 kDa
Shipping	Ice bag