

Anti-Gnao1 Polyclonal Antibody For Wb, If-P, Elisa - P09471

Item Number: CM0023385

Introduction

Rabbit polyclonal antibody against human GNAO1, suitable for WB, IF-P, ELISA. Cross-reacts with mouse and rat. Targets G protein subunit alpha o involved in GPCR signaling and cAMP regulation.

[Learn More](#)

Parameter	Value
Product Name	GNAO1 Rabbit pAb
Remarks/Alias	GNAO; HG1G; DEE17; NEDIM; EIEE17; HLA-DQB1; G-ALPHA-o; GNAO1
Species	Human
GeneID (Human)	2775
GeneID	2775
Immunogen	Recombinant protein corresponding to amino acids 60-180 of human GNAO1 (NP_066268.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09471
Protein Weight	40kDa
Shipping	Ice bag