

Anti-Gnaq Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P50148

Item Number: CM0013332



Introduction

Rabbit polyclonal antibody against human GNAQ for WB, IHC-P, IF/ICC, and ELISA. Detects GNAQ in human, mouse, and rat samples. Suitable for GPCR signaling research.

[Learn More](#)

Parameter	Value
Product Name	GNAQ Rabbit pAb
Remarks/Alias	GAQ; SWS; CMAL; CMC1; G-ALPHA-q; GNAQ
Species	Human
Gene ID (Human)	2776
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-359 of human GNAQ (NP_002063.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50148
Protein Weight	42kDa
Shipping	Ice bag