

Anti-Gnas Rabbit Polyclonal Antibody For Wb, Elisa - O95467 / P63092 / P84996 / Q5Jwf2

Item Number: CM0028960

Introduction

High-quality rabbit polyclonal antibody against human GNAS, validated for Western blot and ELISA. Detects multiple GNAS isoforms including Gs-alpha and NESP55, with cross-reactivity to mouse. Ideal for signal transduction and neuroendocrine research.

[Learn More](#)

Parameter	Value
Product Name	GNAS Rabbit pAb
Remarks/Alias	AHO; GSA; GSP; POH; GPSA; NESP; SCG6; SgVI; GNAS1; PITA3; C20orf45; GNAS
Species	Human
Gene ID (Human)	2778
Gene ID	2778
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GNAS (P63092).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95467/P63092/P84996/Q5JWF2
Protein Weight	28kDa/44kDa/45kDa/77kDa/109kDa/111kDa
Shipping	Ice bag