


Anti-Gnat1 Rabbit Polyclonal Antibody For Wb, Elisa - P11488

Item Number: CM0029637



Introduction

Rabbit polyclonal antibody to GNAT1 (Transducin alpha-1 chain), validated for WB and ELISA, cross-reacts with human, mouse, rat. Suitable for phototransduction and retinal research.

[Learn More](#)

Parameter	Value
Product Name	GNAT1 Rabbit pAb
Remarks/Alias	GBT1; HG1F; GNATR; CSNB1G; CSNBAD3; GNAT1
Species	Human
Gene ID (Human)	2779
Gene ID	2779
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 58-178 of human GNAT1 (NP_000163.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11488
Protein Weight	40kDa
Shipping	Ice bag